# SHIATSU- MASSAGEKISSEN

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## Guide

nah  vital  
just. care.

**DE**

ab Seite 4

**EN**

from page 13

**FR**

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**ES**

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**IT**

da pagina 39

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Output: 12 V DC

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## Dear Customer,

Relaxation benefits the body and the soul. Every massage technique attempts to produce a feelgood moment in its own way. Shiatsu massage, a Japanese method that has existed for 100 years, attempts to do this primarily through touches that reach deep into the tissue.

We are delighted that you have chosen the **nah-vital shiatsu massage cushion**. It mimics the shiatsu massage with 4 massage balls that move over the skin in a circular motion. Thanks to its ergonomic shape, it can be used on the neck, the shoulder blades or even the lower back. You can even use the enclosed foot loops, which are easy to attach with a hook-and-loop fastener, to enjoy a pleasant foot massage. With three speed levels and a heat function that can be activated, you can indulge in your own individual feelgood moment – at home or even using the adapter in a car. If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website: **[www.nah-vital.de](http://www.nah-vital.de)**

We hope you have a lot of fun with your **nah-vital shiatsu massage cushion**.

## Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

---

## LIST OF CONTENTS

Items Supplied	14
At A Glance	14
Symbols	15
Signal Words	15
Intended Use	15
When May the Device Not Be Used?	16
Safety Notices	16
Use	18
Cleaning and Storage	19
Troubleshooting	19
Disposal	20
Technical Data	20

## ITEMS SUPPLIED

- 1x massage device
- 1x foot loops
- 1x mains adapter
- 1x car adapter
- 1x operating instructions

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device.

**Never remove the rating plate and any possible warnings!**

## AT A GLANCE

- |  |  |
|--|--|
| <b>1 direction</b> Select constant direction of rotation of the massage balls                              | <b>6</b> Hook-and-loop fastener for the foot loops   |
| <b>2 power</b> Switch the massage balls with automatic change of direction and red light / heat on and off | <b>7</b> Detachable foot loops                       |
| <b>3</b> Elastic fastening strap on the back   | <b>8</b> Mains adapter                               |
| <b>4 heat</b> Switch red light / heat off and on   | <b>9</b> Connecting socket                           |
| <b>5 speed</b> Set rotational speed at three levels  | <b>10</b> 2 x 4 massage balls with red light element |

## SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Alternating current

Direct current

Only use indoors!

Protection class II

Polarity of the barrel connector

Switching power supply

Safe electrical isolation

## SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

**DANGER** – high risk, results in serious injury or death

**WARNING** – medium risk, may result in serious injury or death

**NOTICE** – may result in risk of damage to material

## INTENDED USE

- This device is intended to be used for loosening the neck and back muscles and delivers relaxation. The device may also be used to massage the abdominal, buttock and leg muscles as well as the feet.
- The device is not designed for medical or therapeutic purposes. It is not a substitute for medical treatment. If complaints persist, consult a doctor.
- The device is for personal use only and is not intended for commercial applications.
- Use the device only for the specified purpose and as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- Do not use the device while driving a vehicle.
- Do not use the device if you have a fever or physical illness/ dental problems.
- Do not use the device while you are connected to medical equipment, e.g. an ECG machine.
- Children, people who are sleeping or unconscious and animals must not be treated with the device.
- Do not use the device while driving a vehicle.
- The following people may use the device **only following consultation with a doctor**:
  - people who are fitted with a pacemaker or similar medical implants,
  - people who suffer from osteoporosis,
  - people with heart disease or who suffer from epilepsy,
  - people who suffer from circulation problems or thrombosis,
  - people who have undergone an operation for which the treatment is not yet complete,
  - people under the influence of medication,
  - people who belong to a risk group or who have concerns about their health.
- Do not use the device in the vicinity of or on:
  - implants,
  - the rib cage, close to your heart,
  - regions which have recently been operated on or are scarred,
  - injuries and areas of swelling,
  - acne, birthmarks and warts,
  - injured, inflamed, (sun-)burnt, irritated or very sensitive skin,
  - body parts with no feeling.

If you are unsure whether the device may be used, consult a responsible doctor.

## SAFETY NOTICES

- This device is not intended for use by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely.
- Children must be supervised to ensure that they do not play with the device.
- The device has a heating function. People who are insensitive to heat must be careful when using the heat function.
- Keep sharp and pointed objects away from the device to ensure that it is not damaged.
- The device must always be disconnected from the mains power when it is left unattended and before cleaning it.
- In the interest of your own safety, inspect the device for damage each time before you use it. Only use it if it is not damaged!

- Do not use the mains adapter (CP1225) with another device and never use the device with another mains adapter or car adapter.
- Do not make any modifications to the device or the connecting cable. The device may only be taken apart and / or repaired by the manufacturer, customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards. If the **connecting cable** of this device is damaged, it must be replaced with a special connecting cable which is available from the manufacturer or its customer service department.

### **DANGER – Danger of Electric Shock!**

- Use and store the device only in closed rooms.
- Do not use and store it in rooms with high humidity.
- Never immerse the device and the connecting cable in water or other liquids and ensure that these components cannot fall into water or become wet.
- If the device falls into water while it is connected to the mains power, pull the plug out of the plug socket immediately! **Only after this** should you retrieve the device from the water.
- Never touch the device and the connecting cable with wet hands when the device is connected to the mains power.
- Do not use the device if it has malfunctioned or fallen into water. Have it checked in a specialist workshop before using it again.

### **WARNING – Danger of Injury!**

- Keep the packaging material away from children and animals. There is a danger of suffocation.
- Make sure that the connecting cable is always kept out of the reach of children and animals. There is a danger of strangulation.

- The device is **not** a toy! Do not leave children and animals unsupervised with the device. There is a danger of injury from incorrect use.
- Do not insert any fingers or other body parts into the device.
- Make sure that the connecting cable does not present a trip hazard when it is connected.

### **WARNING – Danger of Fire**

- Do not use the device under a blanket or similar item or in the vicinity of combustible material.

### **NOTICE – Risk of Damage to Material and Property**

- Use and store the device at room temperature. Do not expose it to extreme temperatures, temperature fluctuations, fire, heat sources and direct sunlight for a prolonged period.
- Only connect the mains adapter to a plug socket that is properly installed and matches the rating plate of the mains adapter. The plug socket must also be readily accessible after connection so that the connection to the mains can quickly be isolated.
- Only use suitable extension cables.
- Protect the device from impacts.

- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- In the event of damage / faults (e.g. unusual noises, odour or smoke development) and in a thunderstorm, switch the device off immediately and disconnect it from the power source.
- Never pull the mains cable to disconnect the car adapter or mains adapter from the mains power. Never pull or carry the device by the connecting cable.
- Make sure that the connecting cable is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- The device must be switched off when inserting or removing the mains plug.

## USE

### Please Note!

- Do not use the device continuously for longer than 30 minutes. Do not massage one group of muscles for longer than 15 minutes.
  - If you experience an unpleasant sensation while you are using the device, stop performing the massage immediately!
  - No pressure should be exerted on the device as this will slow down or stop the massage and the device could be damaged. Do not lie or stand on the device.
  - For reasons of hygiene, avoid direct contact of the skin with the device wherever possible. It should not come into contact with skin that has had cream or oil applied to it.
- A tingling sensation / itching during or after the massage is normal due to the circulation-boosting effect.
  - The device should not be used before going to sleep. The massage has a stimulating effect, which may lead to problems getting to sleep.

We recommend that you wait an hour after mealtimes before using the device.

When this device performs a massage, four massage heads massage the muscles twice. This is intended to relax them and boost circulation. The device is designed in particular for use on the upper back area, in particular the neck. However, it may also be used to massage abdominal, buttock and leg muscles. Furthermore, the device also massages the soles of the feet. To do this, the foot loops that are supplied with the device (7) are attached to the hook-and-loop fasteners (6).

1. First plug the mains adapter (8) or the car adapter into the connecting socket (9) and then into an easily accessible plug socket or the cigarette lighter in a vehicle. When it is used in a vehicle, turn the ignition key to supply power to the device.
2. If the feet are to be massaged, fix the foot loops (7) to the hook-and-loop fasteners (6) and, while sitting down, place your feet inside them.
3. Otherwise place the massage balls (10) onto the part of the body which is to be massaged. For a neck massage, the fastening strap (3) can be placed over a headrest or the back of a chair. It has a hook-and-loop fastener for this purpose.
4. To start the massage, press the **power** button (2). The massage balls rotate and change

direction roughly every minute. The red light is switched on.

5. The massage can be adapted with the other three buttons:
  - set just one direction of rotation or change it with the **direction** button (1)
  - switch the red light / heat on / off with the **heat** button (4)
  - change the rotational speed of the massage balls to three levels with the **speed** button (5)
6. The device switches off automatically after 15 minutes. To switch it off sooner, press the **power** button.
7. Disconnect the connecting cable from the power source and the device.

## CLEANING AND STORAGE

### Please Note!

- Do not use any corrosive or abrasive cleaning agents to clean the device. These may damage the surface.
- 
- Wipe down the device and the foot loops (7) with a damp cloth or sponge and allow them to dry in the air. **Never place on a radiator or dry with a hair dryer!**
  - Do not wind the mains cable around the car adapter or mains adapter.
  - Store the device and accessories in a cool, dry place which is protected from direct sunlight and dust.
  - Store out of the reach of children and animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

### The device does not work or switches itself off prematurely.

The device, car adapter or mains adapter and power source are not connected together correctly.

- ▶ Check and correct all plug connections.

The power source is not supplying any power.

- ▶ Try another power source (plug socket).

- ▶ In a vehicle: Turn the ignition key.

There is too much weight on the device.

- ▶ Reduce the weight on the device.

The device has been operating for too long and has overheated.

- ▶ Disconnect the device from the power source and allow it to cool down. Always disconnect the device from the mains when it is not in use.

### The device gets too warm.

The device has been operating for too long and has overheated.

- ▶ Disconnect the device from the power source and allow it to cool down. Always disconnect the device from the mains when it is not in use.

### The device makes noises.

Operating noises are a normal consequence of the mechanics and not a defect.

- ▶ If the noises increase, switch off the device for approx. 30 minutes.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of this device as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

## TECHNICAL DATA

Article number: 00616  
ID of operating instructions: Z 00616 M NV V1 0521 md

### Device

Model number: C37-MP064  
Voltage supply: 12 V DC  
Power: max. 24 W  
Timer: 15 minutes

### Mains adapter

Manufactured for DS Produkte GmbH

HRB 8937 HL

Am Heisterbusch 1, 19258 Gallin, Germany

Model number: CP1225  
Input: 180 – 240 V AC 50 / 60 Hz; 0.8 A  
Output: 12 V DC; 2.5 A  
Power: 30 W  
Protection class: II  
Type of protection: IP20  
Ø efficiency in operation: approx. 88.08%  
Power consumption on no-load: 0.08 W  
Efficiency on low load (10 %): 86.83%

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All rights reserved.
Prodotto per DS Produkte GmbH