## INHALT

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## CONTENTS

Items Supplied	12	Attachments and Areas of Application	17
At A Glance	12	Intensity Levels	18
Symbols	13	Use	18
Signal Words	13	Care and Storage	19
Intended Use	14	Troubleshooting	19
When May the Device Not Be Used?	14	Technical Data	20
Safety Notices	15	Disposal	21
Charging Up the Battery	16		

## ITEMS SUPPLIED

- Mini massage gun
- Fork attachment
- Flat attachment
- Cartridge attachment
- Sphere attachment
- Heating attachment
- Cooling attachment
- USB cable (A connector to C connector)
- Storage bag
- Operating instructions

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device. **Never remove the rating plate and any possible warnings!**

## AT A GLANCE

(Picture A)

- 1 Mount for the attachment
- 2 Main unit
- 3 Handle
- 4 Intensity level control lamps
- 5  Start / stop / set intensity level
- 6 Charging socket
- 7 Charge control lamp
- 8 Cartridge attachment
- 9 Flat attachment
- 10 Fork attachment
- 11 Sphere attachment
- 12 Heating attachment
- 13 Cooling attachment

Key numbers are shown as follows:  
(1)

Picture references are shown as follows: (Picture A)

## Dear Customer,

We are delighted that you have chosen our **MAXXMEE Mini Massage Gun Hot & Cold Pro**. It operates at up to 3200 strokes per minute and thus delivers a superb powerful massage. The clever bit is that an integrated sensor changes the number of strokes per minute if you slightly increase the pressure while using the device – so you can experience an even more intense massage. Thanks to the 6 different attachments, you can customise the massage individually to suit your needs. Either choose one of the attachments for specific muscle groups or use the heating or cooling attachment to vigorously stimulate blood flow. If you have any questions about the device and about spare parts / accessories, contact the customer service department via our website:

[www.ds-group.de/kundenservice](http://www.ds-group.de/kundenservice)

We hope you have a lot of fun with your **MAXXMEE Mini Massage Gun Hot & Cold Pro**.

## Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

## SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Direct current

## SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

**WARNING** – medium risk, may result in serious injury or death

**NOTICE** – may result in risk of damage to material

## INTENDED USE

- The device is intended to be used for massaging the whole body and is used for relaxation. **It is not a substitute for medically required treatment.**
- The device is **not** a toy.
- The heating and cooling attachments are not suitable for people who are insensitive to the cold or heat.
- The device should not be used before going to sleep. Massages have a stimulating effect, which may lead to problems getting to sleep.
- The device is for personal use only and is not intended for commercial or medical applications.
- The device should only be used as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- Children, people who are sleeping or unconscious, other vulnerable people who are unable to react to overheating or undercooling and animals must not be treated with the device.
- The following people may use the device **only following consultation with a doctor**:
  - people who are fitted with a pacemaker or similar medical implants;
  - people with heart disease or who suffer from epilepsy;
  - people who have back complaints or spinal disorders;
  - people who suffer from osteomyelitis, bone inflammation (in particular if the spinal column is damaged) or from osteoporosis;
  - women who are pregnant or have just given birth;
  - people who suffer from circulation problems or thrombosis;
  - people who have undergone an operation for which the treatment is not yet complete;
  - people with physical abnormalities;
  - people who have a malignant tumour.
- Do not use the device in the vicinity of or on:
  - injuries,
  - swelling,
  - skin irritations,
  - infected skin,
  - warts,
  - fungal infections,
  - abnormal changes in the skin,
  - bones,
  - the head,
  - the genital area,
  - body parts which have been anaesthetised or have no feeling.

- **Only trained people may use the device on joints or around joints.**

If you are unsure whether the device may be used, consult a responsible doctor.

## SAFETY NOTICES

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and /or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- Do not make any modifications to the device. It may only be taken apart and / or repaired by the manufacturer, customer service department or a similarly qualified person (e. g. specialist workshop) in order to avoid any hazards.
- If the battery leaks, avoid physical contact with the battery acid. Wear protective gloves if necessary. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately. There is a danger of burns from battery acid.
- In the interest of your own safety, inspect the device, the attachments and the USB cable for damage each time before you use it. Use the device only if it is in working order.
- Make sure that the USB cable does not present a trip hazard.
- Do not use the device or the heating attachment if it has malfunctioned or fallen into water. Have it checked in a specialist workshop before using it again.

### WARNING – Danger of Injury

- Keep children and animals away from the packaging material. There is a danger of suffocation.
- Make sure that the USB cable is always kept out of the reach of small children and animals. There is a danger of strangulation.
- **WARNING – Danger from Rechargeable Batteries**
- The main unit and the heating attachment, each with a built-in battery, must not be taken apart, thrown into a fire or exposed to high temperatures of over +40 °C. The battery could ignite or explode.
- Make sure that the device, and in particular the main unit and the heating attachment, cannot become wet. Do not use the device and the heating attachment if liquid has got inside. The battery in the main unit or in the heating attachment could ignite or explode.

- Do not cover the device or the heating attachment during charging to prevent overheating.

### NOTICE – Risk of Damage to Material and Property

- The device or the heating attachment must be switched off when they are connected to or disconnected from the mains power.
- Make sure that the connecting cable is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- When the USB mains adapter is pulled out of the plug socket, always pull on the USB mains adapter and never the USB cable.
- Do not pull or carry the device by the USB cable.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- Protect the device from heat sources, naked flames, sub-zero temperatures, persistent moisture, wet conditions and impacts.

## CHARGING UP THE BATTERY

### Please Note!

- A USB mains adapter corresponding to the technical data of the device or the heating attachment is also required (see rating plate or “Technical Data”). This is not supplied with the device.
- Risk of a short circuit! Do not stick any conducting objects (e.g. a screwdriver or similar implement) into the charging sockets.
- Charge the device and the heating attachment only in closed rooms.
- Do not charge the device or the heating attachment in rooms with high humidity.

It is essential that you charge up the battery of the main unit (2) and the heating attachment (12) fully before you first use the device in order to achieve the maximum charging capacity.

1. Plug the USB-C connector of the USB cable into the charging socket (6) on the device or on the back of the heating attachment and the USB connector into the socket of a suitable USB mains adapter.
2. Connect the USB mains adapter to a plug socket that is easily accessible. The charge control lamp (7) on the device or heating attachment lights up

red. When it lights up blue (main unit) or goes out (heating attachment), the battery is fully charged.

3. Pull the USB mains adapter out of the plug socket and disconnect the USB cable from the USB mains adapter and the device or heating attachment.

## ATTACHMENTS AND AREAS OF APPLICATION

(Pictures C–H)

- **Cartridge attachment (8):**  
for deeper tissue and for the feet, hands etc.
- **Flat attachment (9):**  
for relaxing the whole body
- **Fork attachment (10):**  
for the neck, back extensors, calves or similar muscles
- **Sphere attachment (11):**  
for large muscle groups such as quadriceps and gluteus maximus and for targeted relaxation for sore muscles
- **Heating attachment (12):**  
for treating tired muscles
- **Cooling attachment (13):**  
for treating tense muscles

### Heating and Cooling Attachment

In contrast to the other attachments, the heating and cooling attachment (12+13) must be prepared and can be adjusted.

The **cooling attachment** must be placed in the fridge approx. 30 minutes prior to use.

The **heating attachment** has three temperature levels. The level can be identified from the colour ring on the attachment. To activate the heat, press the button on the back of the heating attachment for approx. 2 seconds. The attachment heats up to the lowest level (40 °C) and the colour ring lights up green. To switch to the next level, briefly press the same button (blue: 45 °C > red: 50 °C).

To deactivate the heat, press the button again for 2 seconds; the colour ring goes out.

The heat is automatically deactivated after 15 minutes.

## INTENSITY LEVELS

The device can be operated on 6 different intensity levels. In addition, it has a pressure sensor, which means as soon as a little pressure is

exerted on the device, the number of strokes increases again within the intensity level that is set.

Level	Normal	With pressure
1	1350 strokes / minute	1800 strokes / minute
2	1700 strokes / minute	2100 strokes / minute
3	1950 strokes / minute	2400 strokes / minute
4	2200 strokes / minute	2700 strokes / minute
5	2450 strokes / minute	3000 strokes / minute
6	2750 strokes / minute	3200 strokes / minute

## USE

### Please Note!

- Do not use the device under a blanket or similar cover to prevent overheating.
- Do not use the device when it is connected to the mains power.
- For hygiene reasons, the sphere attachment should not be used on skin which has had cream or oil applied.
- **If you use the heating attachment, keep moving it and do not hold it in one place for too long.**
- If you experience an unpleasant feeling or tingling sensation in the hand holding the device, stop using it and consult a doctor if necessary.

Before using the device, carefully read the "Safety Notices" and "When May the Device Not Be Used?" chapters and familiarise yourself with how to use the device. **However, reading the operating instructions is no substitute for expert advice from a doctor!**

1. Insert the attachment you want as far as it will go into the mount (1) on the main unit (2).
2. Press the  button (5) for approx. 3 seconds to switch on the device. The charge control lamp (7) briefly lights up red and blue. The control lamps (4) briefly light up one after the other. The device starts to operate on intensity level 1.

- Press the  button briefly to switch between the intensity levels. The number of control lamps lit indicates which intensity level is currently set.

We recommend starting with the lowest level and gradually increasing the level of intensity.

- Hold the device by the handle (3) and guide the attachment over the area of the body you want.

**Please Note:** Even though the intensity can be increased with a small amount of pressure, you should **not exert too much pressure!**

We recommend that you do not massage one area of muscle for longer than max. 15 minutes in one go.

- To switch off the device, press the  button for approx. 3 seconds. The control lamps go out.

## CARE AND STORAGE

### Please Note!

- Do not use any caustic or abrasive cleaning agents or cleaning pads to clean the device. These may damage the surfaces.

**Make sure that the device is switched off and is not connected to the mains power.**

- Remove the used attachment from the main unit (2) and clean it. All attachments **apart from the sphere attachment (11) and the heating attachment (12)** can be cleaned with warm water and a cleaning agent.
- Wipe down the main unit, the sphere attachment and the heating attachment with a slightly damp cloth.
- Store the device and all accessories in the storage bag in a dry place which is protected from sunlight and is out of the reach of children and animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

The charge control lamp (7) of the device or the heating attachment does not light up red when the device or the heating attachment are connected to the mains power.

The USB mains adapter is not plugged correctly in the plug socket or the USB connector is not plugged correctly in the charging socket (6) of the device or the heating attachment or in the connecting socket on the USB mains adapter.

- Correct the fit of the plug connections.

EN

The plug socket is defective.

- ▶ Try another plug socket.

There is no mains voltage.

- ▶ Check the fuse of the mains connection.

The battery is fully charged.

- ▶ Disconnect the device from the mains power.

The device does not start.

The battery is empty.

- ▶ Charge up the battery.

The heating attachment (12) does not work. / The colour ring of the heating attachment flashes and then goes out.

The battery is empty.

- ▶ Charge up the battery.

## TECHNICAL DATA

Article number: 02389  
 ID of operating instructions: Z 02389 M DS V2 0322 dk

### Device

Model number: BX-710  
 Battery: 11.1 V (3x Li-ion battery, 1800 mAh)  
 Voltage supply: 5 V, 1 A USB mains adapter (not supplied with the device)

Battery charge time: approx. 8 hours  
 Battery operating time: up to approx. 5 hours (with fully charged battery)

Strokes per minute: max. 3200

### Heating attachment

Battery: 3.7 V (1 x Li-ion battery, 650 mAh)  
 Voltage supply: 5 V, 1 A USB mains adapter (not supplied with the device)

Temperature levels: green: 40 °C, blue: 45 °C, red: 50 °C  
 Battery charge time: approx. 2 hours  
 Battery operating time: up to approx. 75 minutes (with fully charged battery)

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

(Rechargeable) batteries must not be disposed of along with the domestic waste.

Consumers are legally obliged to dispose of (rechargeable) batteries separately.

(Rechargeable) batteries can be handed in free of charge at a collection point of your local authority / district or in retailers so that they can be disposed of in an environmentally friendly way and valuable raw materials can be recovered. If they are not disposed of properly, toxic ingredients may escape into the environment and cause harm to the health of humans, animals and plants.

(Rechargeable) batteries contained in electrical devices must be disposed of separately from them whenever possible. Only hand in (rechargeable) batteries when they are discharged.

In the case of (rechargeable) batteries containing lithium, tape over the poles before you dispose of them to prevent an external short circuit. A short circuit may cause a fire or explosion.

The adjacent symbol (crossed-out bin underlined) means that used devices should not be disposed of in the domestic waste, but should instead be taken to special collection and return systems.

Owners of used devices must remove from the used device any used batteries or rechargeable batteries which are not encapsulated by the used device and can be removed in a non-destructive way and dispose of them separately (see also the section on battery disposal).

Owners of used devices from private households can take them free of charge to the collection points of the public waste disposal authorities or the disposal points which are set up by manufacturers and distributors under the Electrical Equipment Act. Businesses which market electrical and electronic devices are also obliged to take back used devices.

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All rights reserved.