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## LIST OF CONTENTS

Items Supplied	11	Attachments and Areas of Application	17
At A Glance	11	Intensity Levels	17
Symbols	12	Use	17
Signal Words	12	Care and Storage	18
Intended Use	13	Troubleshooting	18
When May the Device Not Be Used?	13	Disposal	19
Safety Notices	14	Technical Data	20
Charging Up the Battery	16		

## ITEMS SUPPLIED

- Massage gun
- Fork attachment
- Flat attachment
- Cartridge attachment
- Sphere attachment
- USB cable (A-connector to USB C-connector)
- Operating instructions

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device.

**Never remove the rating plate and any possible warnings!**

For charging from the mains, a suitable USB mains adapter is also required (see the "Technical Data" chapter).

## AT A GLANCE

(Picture A)

- 1 Mount for the attachment
- 2 Main unit
- 3 Handle
- 4 Intensity level control lamps
- 5 Start / stop / set intensity level
- 6 Charging socket
- 7 Sphere attachment
- 8 Fork attachment
- 9 Flat attachment
- 10 Cartridge attachment

Key numbers are shown as follows: (1)

Picture references are shown as follows: (Picture A)

## Dear Customer,

We are delighted that you have chosen our VITALmaxx Massage Gun Power Grip. It scores points not only with its intense power boost massage, but also its clever shape which delivers even more grip strength. Depending on the muscle group, even areas which are difficult to reach receive an ideal massage and workout. The compact massage gun operates at up to 3200 strokes per minute and thus delivers a superb powerful massage. If you have any questions about the device and about spare parts / accessories, contact the customer service department via our website:

[www.ds-group.de/kundenservice](http://www.ds-group.de/kundenservice)

## Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

## SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Direct current

## SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

**WARNING** – medium risk, may result in serious injury or death

**NOTICE** – may result in risk of damage to material

## INTENDED USE

- The device is intended to be used for massaging the whole body and is used for relaxation. **It is not a substitute for medically required treatment.**
- The device is **not** a toy.
- The device should not be used before going to sleep. Massages have a stimulating effect, which may lead to problems getting to sleep.
- The device is for personal use only and is not intended for commercial or medical applications.
- The device should only be used as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- Children, people who are sleeping or unconscious and animals must not be treated with the device.
- The following people may use the device **only following consultation with a doctor:**
  - people who are fitted with a pacemaker or similar medical implants;
  - people with heart disease or who suffer from epilepsy;
  - people who have back complaints or spinal disorders;
  - people who suffer from osteomyelitis, bone inflammation (in particular if the spinal column is damaged) or from osteoporosis;
  - women who are pregnant or have just given birth;
  - people who suffer from circulation problems or thrombosis;
  - people who have undergone an operation for which the treatment is not yet complete;
  - people with physical abnormalities;
  - people who have a malignant tumour.

## EN

- Do not use the device in the vicinity of or on:
  - injuries,
  - swelling,
  - skin irritations,
  - infected skin,
  - warts,
  - fungal infections,
  - abnormal changes in the skin,
  - bones,
  - the head,
  - the genital area,
  - body parts which have been anaesthetised or have no feeling.
- **Only trained people may use the device on or around joints.**

If you are unsure whether the device may be used, consult a responsible doctor.

## SAFETY NOTICES

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- **Children** must not play with the device.
- Cleaning and **user maintenance** must not be carried out by children, unless they are supervised.
- The device must always be disconnected from the mains power when left unattended and before cleaning it.
- Do not make any modifications to the device. If the device or USB cable are damaged, they must be replaced or repaired by the manufacturer, customer service department or a similarly qualified person (e. g. specialist workshop) in order to avoid any hazards.

### DANGER – Danger of Electric Shock

- Charge and store the device only in closed rooms.
- Do not use, charge or store it in rooms with high humidity.
- Never immerse the device and the USB cable in water or other liquids and ensure that they cannot fall into water or become wet.
- If the device falls into water while it is connected to the mains power, switch off the power supply immediately. **Only after this** should you pull the device out of the water!
- Never touch the device and the USB cable with wet hands when these components are connected to the mains power.
- Do not use the device if it has malfunctioned, fallen or been dropped into water. Have it checked in a specialist workshop before using it again.

### WARNING – Danger of Injury

- Keep the packaging material away from children and animals. There is a danger of suffocation.
- Make sure that the USB cable is always kept out of the reach of small children and animals. There is a danger of strangulation.
- If the battery leaks, do not allow your skin, eyes and mucous membranes to come into contact with the battery

acid. Wear protective gloves if necessary. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately. There is a danger of burns from battery acid.

- In the interest of your own safety, inspect the device and the USB cable for damage each time before you use it. Use the device only if it is in working order.
- Make sure that the USB cable does not present a trip hazard.

### WARNING – Danger from Rechargeable Batteries

- The main unit with a built-in battery must not be taken apart, thrown into a fire or exposed to high temperatures of over +40°C. The battery could ignite or explode. Do not store the main unit in a vehicle where it might be exposed to extreme temperatures.
- Do not use the device if liquid has got inside it. The battery in the main unit could ignite or explode.
- Do not use the device if the battery is leaking.

### WARNING – Danger of Fire

- Do not charge up the device in the vicinity of highly flammable materials (curtains, textiles, etc.).
- Do not cover the device during operation or charging in order to prevent it from overheating.

## NOTICE – Risk of Damage to Material and Property

- Make sure that the USB cable is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- Fully unwind the USB cable before connecting a USB mains adapter.
- When the USB mains adapter is pulled out of the plug socket, always pull on the USB mains adapter and never the USB cable.
- Do not pull or carry the device by the USB cable.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- Do not place any heavy objects on the device or the USB cable.
- Protect the device from heat sources, flames, sub-zero temperatures, persistent moisture and impacts.

## CHARGING UP THE BATTERY

### Please Note!

- Only connect the USB mains adapter to a plug socket that is properly installed and matches its technical data. The plug socket must be readily accessible after connection so that the connection to the mains can quickly be isolated.
- Do not charge up on a computer because the computer might be damaged by the high charging current.
- Pull the mains adapter out of the socket if an error occurs during charging or before a thunderstorm.

It is essential that you charge up the battery fully before you first use the device in order to achieve the maximum charging capacity.

The control lamps (4) flash when the battery charge status drops below 5%.

1. Plug the USB-C connector of the USB cable into the charging socket (6) on the device and the USB-A connector into the socket of a suitable USB mains adapter.
2. Connect the USB mains adapter to a plug socket. The control lamps flash. When they flash constantly, the battery is fully charged.
3. Pull the USB mains adapter out of the plug socket and disconnect the USB cable from the USB mains adapter and the device.

## ATTACHMENTS AND AREAS OF APPLICATION

(Pictures C–F)

- **Sphere attachment (7):**  
for large muscle groups such as quadriceps and gluteus maximus and for targeted relaxation for sore muscles
- **Fork attachment (8):**  
for the neck, back extensors, calves or similar muscles
- **Flat attachment (9):**  
for relaxing the whole body
- **Cartridge attachment (10):**  
for deeper tissue and for the feet, hands etc.

## INTENSITY LEVELS

The device can be operated on 5 different intensity levels. In addition, it has a pressure sensor, which means as soon as a little pressure is exerted on the device, the number of strokes increases again within the intensity level that is set.

Level	without pressure
1	1350 strokes / minute
2	1450 strokes / minute
3	1550 strokes / minute
4	1650 strokes / minute
5	1750 strokes / minute

## USE

### Please Note!

- Do not use the device when it is connected to the mains power.
- For hygiene reasons, the sphere attachment should not be used on skin which has had cream or oil applied.
- If you experience an unpleasant feeling or tingling sensation in the hand holding the device, stop using it and consult a doctor if necessary.

Before using the device, carefully read the "Safety Notices" and "When May the Device Not Be Used?" chapters and familiarise yourself with how to use the device. **However, reading the operating instructions is no substitute for expert advice from a doctor!**

1. Insert the attachment you want as far as it will go into the mount (1) on the main unit (2).
2. Press the  button (5) for approx. 3 seconds to switch on the device. The device starts to operate on intensity level 1.
3. Press the  button briefly to switch between the intensity levels. The number of control lamps lit indicates which intensity level is currently set.

We recommend starting with the lowest level and gradually increasing the level of intensity.

4. Hold the device by the handle (3) and guide the attachment over the area of the body you want.

**Please Note:** Even though the intensity can be increased with a small amount of pressure, you should not exert too much pressure!

We recommend that you do not massage one area of muscle for longer than max. 15 minutes in one go.

5. To switch off the device, press the  button for approx. 3 seconds. The control lamps go out.

## CARE AND STORAGE

### Please Note!

- Do not use any sharp or abrasive cleaning agents or cleaning pads to clean the device. These may damage the surfaces.

**Make sure that the device is switched off and is not connected to the mains power.**

- Remove the used attachment from the main unit (2) and clean it. All attachments **apart from the sphere attachment (7)** can be cleaned with warm water and a cleaning agent.
- Wipe down the main unit and the sphere attachment with a slightly damp cloth.

- Store the device and all accessories in a dry place which is protected from sunlight and is out of the reach of children and animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

The control lamps (4) do not light up when the device is connected to the mains power.

The USB mains adapter is not plugged correctly in the plug socket or the USB connector is not plugged correctly in the charging socket (6) or in the connecting socket on the USB mains adapter.

- ▶ Correct the fit of the plug connections.

The plug socket is defective.

- ▶ Try another plug socket.

There is no mains voltage.

- ▶ Check the fuse of the mains connection.

**The device does not start.**

The battery is empty.

- ▶ Charge up the battery.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

(Rechargeable) batteries must not be disposed of along with the domestic waste.

Consumers are legally obliged to dispose of (rechargeable) batteries separately.

(Rechargeable) batteries can be handed in free of charge at a collection point of your local authority/district or in retailers so that they can be disposed of in an environmentally friendly way and valuable raw materials can be recovered. If they are not disposed of properly, toxic ingredients may escape into the environment and cause harm to the health of humans, animals and plants.

(Rechargeable) batteries contained in electrical devices must be disposed of separately from them whenever possible. Only hand in (rechargeable) batteries when they are discharged. In the case of (rechargeable) batteries containing lithium, tape over the poles before you dispose of them to prevent an external short circuit. A short circuit may cause a fire or explosion.

The adjacent symbol (crossed-out bin underlined) means that used devices should not be disposed of in the domestic waste, but should instead be taken to special collection and return systems.

Owners of used devices must remove from the used device any used batteries or rechargeable batteries which are not encapsulated by the used device and can be removed in a non-destructive way and dispose of them separately (see also the section on battery disposal).

Owners of used devices from private households can take them free of charge to the collection points of the public waste disposal authorities or the disposal points which are set up by manufacturers and distributors under the Electrical Equipment Act so that they can be disposed of in an environmentally friendly way and valuable raw materials can be recovered. If they are not disposed of properly, toxic ingredients may escape into the environment and cause harm to the health of humans, animals and plants. Businesses which market electrical and electronic devices are also obliged to take back used devices.

Waste prevention makes an even more valuable contribution to protecting the environment. This means if possible, as well as continuing to use or repair a device yourself, passing it on to another user also represents an ecologically beneficial alternative to disposal.

## TECHNICAL DATA

Article number:	03764 (black); 11039 (white)
ID of operating instructions:	Z 03764_11039 M DS V1.3 0723 as
Model number:	D72-MG003
Battery:	7.4 V DC (1 x Li-ion battery, 1800 mAh)
Power:	25 W
Power supply:	5 V, 1 A USB mains adapter (not supplied with the device)
Battery charge time:	approx. 4 hours
Battery operating time:	approx. 4 to 5.5 hours (with fully charged battery at the lowest intensity level)
Strokes per minute:	max. 3000

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All rights reserved.