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LIST OF CONTENTS

Items Supplied	9
Symbols	10
Signal Words	10
At A Glance	10
Intended Use	11
Safety Instructions	11
Use	14
Cleaning and Storage	14
Troubleshooting	15
Disposal	15
Technical Data	16

Dear Customer,

We are delighted that you have chosen the **massage seat cushion**. If you have any questions about the device and about spare parts / accessories, contact the customer service department via our website: www.ds-group.de/kundenservice

Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

ITEMS SUPPLIED

- 1 x massage seat with controller
- 1 x mains adapter
- 1 x operating instructions

Check the device and the connection for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device. **Never** remove the rating plate and any possible warnings!

SYMBOLS

Danger symbols: These symbols indicate dangers of injury. Read through the associated safety notices carefully and follow them.

Supplementary information

Read operating instructions before use.

Protection class II

Protection class III

Polarity of the coaxial connector

Do not pierce the device with needles.

Do not wash.

Only use indoors!

Direct current

Energy efficiency class V

SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

DANGER – high risk, results in serious injury or death

WARNING – medium risk, may have serious injury or death

CAUTION – low risk, may result in minor or moderate injury

NOTICE – may result in risk of damage to material

AT A GLANCE

Picture A

- 1 Upper back massage zone
- 2 Lower back massage zone
- 3 Hips massage zone with heat that can be activated
- 4 Thighs massage zone
- 5 Storage bag for the controller (not shown)
- 6 Connecting socket for the mains adapter cable (not shown)
- 7 Controller

Picture B

- 8 **WÄRME (HEAT)** Switch the heat on and off
- 9 **INT** Set the massage strength. There are 4 intensity levels.
- 10 **MODUS (MODE)** Switch forward to the next massage programme. There are a total of 5 massage programmes, including **AUTO**, in which the 4 different types of massage are run through automatically.
- 11 **OBEN (UPPER) / UNTEN (LOWER) / HÜFTE (HIPS) / BEIN (LEG)** Deactivate and reactivate individual massage zones.
- 12 **TIMER** Select the duration of the massage. After 15, 30 or 60 minutes, the device switches off automatically.
- 13 **LEDs** indicate the activated massage zones.
- 14 **PWR** Switch the device on and off. After switch-on, the **AUTO** programme is automatically activated for all massage zones.

INTENDED USE

- The device is intended for massaging and for warming your back / buttocks and thighs and is used for relaxation. It is not a replacement for medically required treatment.
- Never stand on the device.
- The device is for personal use only and is not intended for commercial applications.
- The device should only be used as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

SAFETY INSTRUCTIONS

WARNING: Note all safety notices, instructions, illustrations and technical data provided with this device.

Failure to follow the safety notices and instructions may result in electric shock, fire and/or injuries.

- **Children** younger than three years of age must not use this device as they are not capable of reacting to overheating.
- The device must not be used by young **children** aged three or over unless the regulating and / or control units have been preset by a parent or a supervisor or the child has received sufficient instruction on how to operate the regulating and / or control units safely.
- This device may be used by **children** aged 8 and over and by people with reduced physical, sensory or mental abilities or a lack of experience and/or knowledge if they are supervised or have received instruction on how to use the device safely and have understood the dangers resulting from failure to comply with the relevant safety precautions.
- Children must **not** play with the device.
- Cleaning and **user maintenance** may **not** be performed by **children**, unless they are older than 8 years of age and are supervised.
- When the heat function is switched on, the device has a hot surface. People who are not sensitive to heat must be careful when using it.

- The device may not be used by vulnerable people who are unable to react to overheating.
- The device is not intended for use in hospitals.
- Do not use the device if it is wet or damp. It must be completely dry before you start using it.
- No safety pins or other sharp or pointed objects must be attached to the device or inserted into it.
- The device should not be switched on in a folded condition. It must always be spread out prior to use to prevent the build-up of heat.
- Always check the device and accessories for damage or wear before putting them into operation. It may only be used if it is undamaged and does not display any signs of wear and tear.
- Do not make any modifications to the device.
- If the device or the connection is damaged, they must be replaced by the manufacturer, the customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards.
- The device must always be disconnected from the mains power when left unattended and before cleaning it.
- The device should always be allowed to cool down completely before it is washed, folded up and/or stored away.
- To ensure that the device is not damaged, fold up it loosely, do not press it together firmly or bend it. Do not place any objects, animals or people on the folded device.
- Only ever use the device with the mains adapter (XH1200-0800wG) that is supplied with it.

DANGER – Danger of Electric Shock

- Use and store the device only in closed rooms.
- Do not operate or store the device in rooms with high humidity.
- Never immerse the device in water or other liquids! Ensure that the device and the connection cannot fall into water or become wet.
- If the device becomes wet or falls into water while it is connected to the mains power, pull the mains plug out of the plug socket immediately. **Only after this** should you pull the device out of the water. Never attempt to pull it out of the water while it is still connected to the mains power!
- Never device and the connection with wet hands when the device is connected to the mains power.

- Do not use the device if it has malfunctioned or has been dropped in water in order to avoid any hazards. Have it checked in a specialist workshop before using it again.

WARNING – Danger of Fire

- Do not operate the device in the vicinity of combustible material.
- Do not place any combustible materials or materials that melt easily on the device.

WARNING – Danger of Injury

- Keep packaging material well away from children and animals. There is a danger of suffocation.
- Make sure that the connection is always kept out of the reach of small children and animals. There is a danger of strangulation.
- People who have no feeling in their neck or back or people with circulation problems may not use the device.
- Babies and sleeping or unconscious people may **not** be placed on the device.
- **Only use the device if you have no physical illnesses or ailments.** If you are unsure whether the device may be used, consult a doctor.
- People with a pacemaker may use the device only after consultation with their doctor.
- Pregnant women or those who have recently given birth should not use the device, or should use it only after consultation with their doctor.
- Do not use the device near to or on injuries, swelling, skin irritations, infected skin, warts, fungal infections, abnormal changes in the skin or body parts with no feeling.

- Switch off the device immediately if you experience pain or feel unwell during use.

NOTICE – Risk of Damage to Material and Property

- Make sure that the connection is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- Disconnect the device from the mains power if a fault occurs or during a thunderstorm.
- When you disconnect the device from the mains power, always pull on the mains plug and not the cable. Never pull or carry the device by the connection.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- Use and store the device at room temperature, do not expose it to extreme temperatures or temperature fluctuations.
- Protect the device from impacts.
- Do not use the connection with a different device.

USE

Please note!

- The mains plug should only be connected to a plug socket that has been properly installed. The socket must also be readily accessible after connection. The mains voltage must match the technical data for the device (see the "Technical Data" chapter).
- Make sure that the connection does not present a trip hazard.
- Unwind the mains cable fully before connecting the mains adapter.
- To prevent overheating, do not use the device under a blanket or similar cover.
- Do not use it in bed. Do not fall asleep during use.
- The device must not be folded up or covered when it is switched on. Do not place the connection or the controller on or beneath the device during operation.

Do not use the device before going to bed. The massage has a stimulating effect, which may lead to problems getting to sleep.

1. Place the massage seat cushion on a chair, for example, and pull the fastening elastic bands over the backrest so that the seat cushion does not slip.
2. Connect the mains cable to the connecting socket on the device and insert the mains plug of the adapter into an easily accessible plug socket.
3. Sit down on the massage seat cushion and lean back.

4. Take the controller in your hand and press the **PWR** button. The **AUTO** massage programme will then start. By pressing the buttons on the controller, you can set the massage to your liking and activate the heat function (see the "Device Overview" chapter – "Controller"). You can also switch on the heat function exclusively by pressing the **WÄRME (HEAT)** button instead of the **PWR** button.
5. Set the massage / heating duration (15, 30 or 60 minutes) with the **TIMER** button. We recommend that you choose 15 minutes.
6. You can keep the controller inside the pocket on the side of the device during use.
7. When the set time has elapsed, the device switches off automatically. You can also switch it off at any time by pressing the **PWR** button.
8. Remove the mains adapter from the plug socket and disconnect the cable from the device.
9. Allow the device to cool down if the heat function was activated before you put it away or clean it (see "Cleaning and Storage" chapter).

CLEANING AND STORAGE

Please note!

- Do not use any sharp or abrasive cleaning agents or cleaning pads to clean the device. These could damage the surfaces.
-
- Make sure that the device is not connected to the mains power and has cooled down.
 - Wipe down the massage seat cushion and the controller with a damp cloth. Then dry it with a soft cloth or allow it to dry in the air.

- Loosely fold up the massage seat cushion and place it together with the controller in the associated packaging. Store the device in a cool, dry place which is protected from direct sunlight and is out of the reach of children and animals.
- Do not place any objects on the folded massage mat, so that it is not compressed.

TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

Do not attempt to repair an electrical device yourself!

The device does not work.

The device is not supplied with power.

- ▶ Check whether
 - the barrel connector is inserted correctly in the connecting socket.
 - the mains adapter is inserted in the plug socket.
 - the plug socket is defective.
 - the fuse of the mains connection is inserted.

The **PWR** button was not (correctly) pressed.

- ▶ Press the **PWR** button.

All massage zones are deactivated on the controller.

- ▶ Press the button for a massage zone to activate it. The associated LED ought to light up.

No heat is produced.

The device is not switched on.

- ▶ Press the **PWR** button. The LED must light up.

The heat function is not switched on.

- ▶ Press the **WÄRME (HEAT)** button. The LED must light up.

DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

The adjacent symbol (crossed-out bin underlined) means that used devices should not be disposed of in the domestic waste, but should instead be taken to special collection and return systems.

Owners of used devices from private households can take them free of charge to the collection points of the public waste disposal authorities or the disposal points which are set up by manufacturers and distributors under the Electrical Equipment Act so that they can be disposed of in an environmentally friendly way and valuable raw materials can be recovered. If they are not disposed of properly, toxic ingredients may escape into the environment and cause harm to the health of humans, animals and plants. Businesses which market electrical and electronic devices are also obliged to take back used devices.

Waste prevention makes an even more valuable contribution to protecting the environment. This means if possible, as well as continuing to use or repair a device yourself, passing it on to another user also represents an ecologically beneficial alternative to disposal.

TECHNICAL DATA

Article number: 04877
ID of operating instructions: Z 04877 M DS V11 0624 as

Device

Model number: CF-2602
Voltage supply: 12V DC
Power: 9,6W
Protection class: III

Mains adapter

Manufactured for DS Produkte GmbH
HRB 8937 HL
Am Heisterbusch 1,
19258 Gallin, Germany

Model number: XH1200-0800wG
Input: 100–240 V ~ 50/60 Hz
Output: 12 V DC; 800mA
Protection class: II
Ø efficiency in operation: 74.61 %
Efficiency on low load (25 %): 63.99 %
Power consumption on no-load: 0.15W

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All rights reserved.
Gefabriceerd voor DS Produkte GmbH