Inhalt

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Contents

| | |
|--|----|
| Intended Use _____ | 11 |
| When May the Device Not Be Used? _____ | 11 |
| Safety Notices _____ | 12 |
| Charging Up the Battery _____ | 13 |
| Safety Devices _____ | 14 |
| Functions _____ | 14 |
| Use _____ | 14 |
| Cleaning and Storage _____ | 16 |
| Troubleshooting _____ | 16 |
| Disposal _____ | 17 |
| Technical Data _____ | 17 |

Items Supplied

(Picture A)

- 1 Abdominal muscle belt 1 x
- 2 Extension strap 1 x
- 3 USB cable (A connector to Micro B connector) 1 x
- 4 Operating element incl. battery 1 x
- 5 Remote control incl. battery 1 x
 - Storage bag 1 x
 - Operating instructions 1 x

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any films, stickers or transport protection. **Never remove the rating plate and any possible warnings!**

Key numbers are shown as follows: (1)

Picture references are shown as follows:

(Picture A)

Displays are shown as follows: **M1**

Button designations are shown as follows:
ON /OFF

A USB mains adapter is also required.

Symbols

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Symbol for direct current

Protection class III

Signal Words

Signal words indicate hazards if the associated notices are not followed.

DANGER – high risk, results in serious injury or death

WARNING – medium risk, may result in serious injury or death

NOTICE – may result in risk of damage to material

Dear Customer,

We are delighted that you have chosen the **FitEngine I-50 EMS Abdominal Muscle Belt**.

If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website: www.fit-engine.de

We hope you have a lot of fun with your **FitEngine I-50 EMS Abdominal Muscle Belt**.

Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

Intended Use

- This device is intended to be used for stimulating the skin and the tissue beneath and for electrical muscle stimulation.
 - The device is **not** designed for medical or therapeutic purposes.
 - The device is intended for personal domestic use only and is not intended for commercial applications.
 - The device should only be used as described in the operating instructions. Any other use is deemed to be improper.
 - The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.
- The following people may use the device only following consultation with a doctor:
 - people who are fitted with a pacemaker or similar medical implants;
 - people with heart disease or who suffer from epilepsy;
 - people who have spinal disorders or back complaints,
 - women who are pregnant or have just given birth;
 - women during their period,
 - people who have undergone an operation for which the treatment is not yet complete,
 - people with physical abnormalities,
 - people with a malignant tumour,
 - people under the influence of medication,
 - people who suffer from circulation problems or thrombosis;
 - people with nerve damage or paralysis;
 - people who belong to a risk group or who have concerns about their health.

When May the Device Not Be Used?

- Infants, people who are sleeping or unconscious and animals must not be treated with the device.
- Do not use the device if you have a fever.
- Do not use the device if you have calf pain.
- Do not use the device on:
 - body parts with no feeling;
 - regions which have recently been operated on or are scarred;

- injured, infected, inflamed, (sun-) burnt or very sensitive skin;
- injuries and areas of swelling;
- acne, birthmarks, warts, rosacea;

- Do not use the device on your head, neck, face, genital area or on the rib cage close to your heart.
- If you are unsure whether the device may be used, consult a responsible doctor.

Safety Notices

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- Do not make any modifications to the device. It may only be taken apart and / or repaired by the manufacturer, customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards.

WARNING: Do not use this device near a bathtub, shower or any vessel containing water. Never use while bathing or showering!

Danger – Danger of Electric Shock

- If the operating element falls into water while it is connected to the mains power, pull the USB mains adapter out of the plug socket immediately! **Only after this** should you pull the device out of the water and dispose of it.
- Never touch the device and the connecting cable (USB cable and USB mains adapter) with wet hands when the operating element is connected to the mains power.

WARNING – Danger from Rechargeable Batteries

- Make sure that the operating element and the remote control cannot become wet. **Do not immerse in water!** Do not use the operating element and the remote control if liquid has got inside.
- Do not use the device if it has malfunctioned or fallen into water. Have it checked in a specialist workshop before using it again.
- The operating element and the remote control have a built-in battery. They must not be thrown into a fire or exposed to high temperatures of over +40°C. Do not take batteries apart!
- If the battery leaks, do not allow your skin, eyes and mucous membranes to come into contact with the battery acid. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately.

WARNING – Danger of Fire

- Keep the device away from combustible materials and substances.
- Do not cover the device during operation.

WARNING – Danger of Injury

- Keep children and animals away from the packaging material. There is a danger of suffocation.
- Make sure that the USB cable is always kept out of the reach of young

children and animals. There is a danger of strangulation.

- Children must be supervised to ensure that they do not play with the device. There is a danger of injury from incorrect use.
- Do not use the device if you are attached to high-frequency surgical devices or are in the vicinity of such devices or high-frequency commercial equipment.
- Do not use the device in the immediate vicinity of devices which operate using short waves or microwaves. Maintain a minimum gap of 1 metre.
- The device must not be used near a mobile phone. Do not use a mobile phone while using the device.
- Do not use the device while driving a vehicle or while sleeping.
- Batteries can be deadly if swallowed. Batteries, the remote control and operating element should be kept out of the reach of children and animals. If a battery is swallowed, seek medical assistance immediately. There is a danger of burns from battery acid!
- If you experience an unpleasant feeling while using the device, stop using it immediately!
- Stop using the device if you experience skin redness or discolouration that lasts for more than 24 hours.

NOTICE – Risk of Damage to Material and Property

- Remove a dead battery from the operating element or the remote control immediately.
- There are magnets on the operating element and on the remote control. These may interfere with sensitive electronic devices or magnetic data storage media. They should therefore be kept away from the operating element and the remote control.

- Store the device only in closed rooms.
- Do not expose the operating element and the remote control to any moisture.
- Use and store the device at room temperature. Do not expose it to extreme temperatures, temperature fluctuations, fire, heat sources and direct sunlight for a prolonged period.
- Protect the device and accessories from impacts.
- Do not place heavy objects on the device or its accessories.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- In the event of damage / faults, switch the device off immediately.
- Do not connect the operating element to a plug socket during a thunderstorm.

Charging Up the Battery

PLEASE NOTE!

- Do not charge up on a computer because the computer might be damaged by the high charging current.

Remote Control

If the operating element (4) no longer responds to the remote control (5), the battery of the remote control is too weak and needs to be charged.

To connect the device to the plug socket, a USB mains adapter (5 V and $\geq 1A$; not supplied with the device) is required.

1. First plug the USB cable (3) into the Micro-USB socket on the remote control and then connect to a power source with USB interface (preferably via a USB mains adapter to the plug socket).

The charging process begins and **L** is shown on the display of the remote control.

- As soon as **F** appears on the display, the battery is charged.
Disconnect the USB cable from the power source and then from the USB socket.

Operating Element

If the operating element (**4**) cannot be switched on, the battery is too weak and needs to be charged.

To connect the device to the plug socket, a USB mains adapter (5V and ≥ 1 A; not supplied with the device) is required.

- First plug the USB cable (**3**) into the Micro-USB socket of the operating element and then connect to a power source with USB interface (preferably via a USB mains adapter to the plug socket).

The charging process begins and the charge control lamp on the operating element lights up red.

- As soon as the charge control lamp goes out, the battery is charged.
Disconnect the USB cable from the power source and then from the USB socket.

Safety Devices

The operating element (**4**) has two built-in safety functions:

- The operating element automatically switches off after 10 seconds if no contact is established between the skin and silicon pads in the belt (**1**).
- The operating element automatically switches off after 15 minutes of use.

Functions

EMS stands for electrical muscle stimulation. The device is used to stimulate the abdominal muscles.

Beginners should use the device for no more than 15 minutes a day to prevent overstimulation!

Advanced users can use the device for up to 30 minutes a day. We recommend no more than 2 units per day, with a rest phase of at least 4 hours between the units.

Modes

M1 – Warm-up and cool down

M2 – Intensive training, high level

M3 – Intensive training, medium level

M4 – Advanced training, medium level

M5 – Advanced training, high level

M6 – Complete training, high level

Each mode has 9 different intensity levels.

Use

Please Note!

- Do not use straight after eating.
- The device should not be used before going to sleep. Massages have a stimulating effect, which may lead to problems getting to sleep.
- For reasons of hygiene, the device should only be used by one person and cleaned after each use (see the “Cleaning and Storage” chapter). Keep the device clean.
- The skin may feel warm during the treatment. This is normal and not a cause for concern.
- Do not bend the silicone pads.
- The silicone pads must not come into contact during use with metal parts, e.g. on belts, wrist watches, jewellery or similar items.

- We explicitly point out that with some people skin irritations or hypersensitivity may occur as a result of the electrical muscle stimulation. In this case, stop using the device immediately and consult a doctor if necessary.

1. Thoroughly clean and then dry the area which is to be trained.
2. Place the belt (1) around you so that the silicone pads in the belt rest directly on the skin over the abdominal muscles (**Picture B**). Fix the belt using the hook-and-loop fastener.
3. Attach the operating element (4) to the belt using the press studs. Both press studs must fit securely for the device to work.

The belt is suitable for a waist circumference of 65 – 110 cm. This circumference can be extended up to 140 cm with the extension strap (2).

Operating the Device

The operating element (4) can also be controlled with the remote control (5).

1. Switch on the operating element: press the **+** button on the operating element. **M1** and **0** is shown on the display and the operating lamp lights up green. When the device is switched on, the intensity of the stimulation is set to level **0** and mode **M1**. In total there are 6 different modes, each with 9 intensity levels (see the “Functions” chapter).

To activate the remote control, press the **<** button.

2. Press the **+** button on the operating element or on the remote control repeatedly to increase the intensity to the level where stimulation still feels pleasant. Use the **-** button on the operating element or on the remote control to reduce the intensity if necessary.

3. To switch mode, briefly press the **M** button on the operating element or on the remote control.

After each change of mode, the intensity is automatically set to level **0** and needs to be manually increased again.

4. To switch off the operating element, press the **-** button on the operating element several times or the **>** button on the remote control until the display switches off and the operating lamp goes out.

The operating element automatically switches off after 15 minutes.

For more information and inspiration, have a look at our website:

www.fit-engine.de.

Here you will find plenty of training videos and lots of tips for your own personal training.

What You Should Note When Training With EMS

- The most important thing first: Training with electrical muscle stimulation (EMS) is the perfect accompaniment to your regular sports regime. However, it is not a substitute for traditional training!
- The EMS only boosts your strength and is not a substitute for endurance training. To ensure that you train your whole body in a sensible way, you should also focus on strengthening your stamina.
- You should build up strength not only using EMS. Targeted traditional strength training combined with EMS will boost the effect.
- Particularly if you want to improve your performance in a particular sport through EMS, you should replicate movements that are typical of this sport during the training. This will help you to train your coordination.

- It is particularly important to drink a lot during EMS training. This is because firstly large amounts of water can collect in the muscles during the training. So to ensure you do not experience any circulation problems, you should always remember to drink enough. And secondly EMS training increases the secretion of the enzyme creatine kinase (CK), which is eliminated via the kidneys. Having sufficient fluid in the body aids the kidneys.
- Include sufficient recovery phases between your training sessions. We recommend that you train no more than once or twice a week.
- EMS training is not designed for you to endure pain in order to produce supposed better results. This can damage your health more than it benefits it. You should therefore never train too much and not at an intensity level that is too high.

Cleaning and Storage

PLEASE NOTE!

- After each use, store the device in a place where children cannot use it.
 - No liquid must get into the operating element. Never immerse the operating element in water!
 - Do not use any caustic or abrasive cleaning agents or pads to clean the device. These may damage the surfaces.
-
- Remove the operating element (4) from the belt (1) after use.
 - Wash the belt in accordance with the washing instructions on the washing label.
 - Clean the operating element with a slightly damp, soft cloth.
 - Store all parts in a dry place which is protected from moisture, dirt and

direct sunlight and is out of the reach of children and animals.

Troubleshooting

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

Do not attempt to repair an electrical device yourself!

The operating element (4) does not work. / The display does not show anything.

The battery is empty.

- ▶ Charge up the battery.

The remote control (5) does not work.

No batteries or empty batteries have been inserted.

- ▶ Insert / replace batteries.

The batteries are inserted incorrectly.

- ▶ Pay attention to the correct polarity when inserting the batteries.

The operating element (4) repeatedly switches off.

There is no contact between the skin and silicone pads in the belt (1).

- ▶ Make sure that the silicone pads fit correctly. There must be skin contact. There must not be any hairs, clothing, jewellery etc. between the skin and silicone pads.

The operating element (4) is not attached to the belt (1) correctly.

- ▶ Check whether the operating element is correctly attached to the press studs.

The battery is too weak.

- ▶ Charge up the battery.

The display shows the mode and intensity, but no pulses for muscle stimulation are emitted.

The silicone pads in the belt (1) are not fitted correctly.

- ▶ Keep changing the position of the silicone pads until the spot with the strongest muscle contraction has been found.

The intensity level is set too low.

- ▶ Increase the intensity level.

The silicone pads are dirty.

- ▶ Clean the silicone pads.

The battery is too weak.

- ▶ Charge up the battery.

Disposal

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

Batteries should be removed from the remote control prior to disposal of the device and disposed of separately from the device at suitable collection points.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of the operating element with a built-in battery as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

Technical Data

Article number: 05022

ID of operating instructions: Z 05022 M FE V2 0121 uh

Remote Control

Model number: QRM-01

Voltage supply: 3.7 V DC (1 x Li-ion battery, 150 mAh)

Operating Element

Model number: QSF098

Voltage supply: 3.7V DC (1 x Li-ion battery, 150 mAh)

Charge time: approx. 1 h

Operating time: approx. 150 min; automatic switch-off after 15 min

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All rights reserved.